


## FITUR

- Pengoperasiaan yang mudah
- Layar sentuh
- Ukuran Layar TFT 3.5" (320 x 240)
- Terdapat socket TF card
- Dapat menampilkan grafik gelombang ECG, Respirasi dan IBP & Temperature (optional)
- Terdapat indikator baterai
- Built-in 7,4V/4000mAh rechargeable battery



## SPEKIFIKASI

Model	MED-S400
Normal Sinus Rhythm	ECG reference Lead II, output lead ACG, Amplitude 0,5-1,0-1,5-2,0 mV
Arhythmia	Default 80 bpm, Atrial, Ventricular, Premature Beats
Pacemaker	Pace amplitude 2,0 mV, Pace width 2,0 ms
ECG Performance	Standart adult ECG, standart pediatric ECG, QRS wave
Respiration	Accuracy $\pm 2\%$ , impedance, baseline impedance
Temperature	0°C, 24°C, 37°C, 40°C
N.W	600 gr
Dimension	45 (H) x 189 (L) x 141 (W)
Baterai	7,4V (Lithium Rechargeable)
Input Power Adaptor	9 - 15VDC 1000mA