

# PROMIST 1

Mini Compressor Nebulizer



Elitech<sup>®</sup>  
TECHNOVISION



ecatalogue



**LKPP**

Lembaga Kebijakan  
Pengadaan Barang/ Jasa Pemerintah

SUKSESKAN GERAKAN



*“Nebuliser untuk  
Inhalasi yang efektif...”*



**CPAKB**  
CARA PEMBUATAN  
ALAT KESEHATAN  
YANG BAIK  
HK.02.04/VI/011/2016  
PT. SINKO PRIMA ALLOY



# PROMIST 1




## Mini Compressor Nebulizer



Adult Cup & Mouthpiece



Pediatric Cup & Mouthpiece



Pediatric Mask



Adult Mask



Carry Bag



Air Filter



KEMENTERIAN KESEHATAN REPUBLIK INDONESIA  
KEMENKES RI  
AKD 20403318003



84778/A/0001/UK/En



ISO9001:2015



ISO13485:2016  
MS 181050 R

### SPECIFICATION

Model	PROMIST-1	
Power Input	AC 220V, 50Hz ,	
Daya Terpakai	< 60W	
Normal Operating	Ambient Temperature	+10°C - +40°C
	Relative humidity	30% RH - 85% RH
	Aliran Kecepatan Udara	80kpa - 105kpa
Transport Operating	Ambient Temperature	-20°C - +70°C
	Relative humidity	≤ 95%
	Athmosphere	80kpa - 105kpa
Noise Level	50 dB (A)	
Compressor Pressure Range	35 Psi - 50 Psi (241 kpa - 345 kpa)	
Operating Pressure Range	9-16 Psi (62 kpa - 110 kpa)	
Flow Rate Range	5-8 L/Min	
Particle Size (MMAD)	3 micron	
Cup Capacity	2-6 ml	
Dimensi	155 x140x96 mm	
Berat	1,2 kg	

hubungi

Diproduksi oleh:

**PT. Sinko Prima Alloy**  
Surabaya – Indonesia

BR-CN1/2016/00