



Elitech[®]

TECHNOVISION

SonoTrax-B



- Lampu Belakang LED
- Layar Display 2"
- 3-Model Fungsi
- Terdapat Speaker
- Termasuk earphone
- Probe dapat diganti (tambahan)
- Otomatis shut-off
- 4-Level Indikator Baterai
- Menggunakan Baterai 2 x AA
- Waktu Pengoperasian ± 6 jam



SUKSESKAN GERAKAN
100%
INDONESIA

1 TAHUN
GARANSI RESMI
Free Servis dan Spare Parts

10
upgrade
MEMORY

Saya **Produksi Indonesia**
made in INDONESIA

Pocket Fetal Doppler

Model	SONOTRAX-B	
Baterai	2 "AA" alkaline,	
Normal Operating	Ambient Temperature	±5°C~40°C
	Relative humidity	≤80%
	Aliran Kecepatan Udara Ruangan	86kpa-109kpa
Transport Operating	Ambient Temperature	-10°C~+55°C
	Relative humidity	≤93%
	Athmosphere	86kpa-106kpa
	Daya keluaran audio	0,5 W
Jangkauan ukur FHR	60-240 BPM	
Akurasi	±2 BPM	
Resolusi pengukuran FHR	1 BPM	
Probe Frekuensi	Frekuensi kerja probe dengan 2.3MHZ; Unit kompatibel dengan 2~8MHZ Probe	
Waktu Operasi	6 jam	
Dimensi	135mm(L)x95mm(W)x35mm(H)	
Berat	225 gr	

e catalogue

LKPP
Lembaga Kebijakan Pengadaan Barang/ Jasa Pemerintah

KEMENTERIAN KESEHATAN REPUBLIK INDONESIA
KEMENKES RI
AKD 21102800003

UKAS MANAGEMENT SYSTEMS
043
84778/A/0001/UK/En

ISO 9001
REGISTERED FIRM
ISO9001:2015

ISO13485:2016
MS 181050 R
globalgroup ISO 13485

Manufacturer :
PT. Sinko Prima Alloy.
Surabaya - Indonesia

CPAKB
CARA PEMBUATAN ALAT KESEHATAN YANG BAIK
HK.02.04/VI/011/2016
PT. SINKO PRIMA ALLOY



2MHz ~ 8MHz